

BioReady Taq pac.

Kit Components

Cat#	BSA09S2	BSA09M2
Components	200units	500units
10 × Reaction Buffer (with 15mM MgCl ₂)	500μl	1000μl
BioReady Taq (5U/μl)	100μl	200μl
2.5mM dNTP Mixture	350μl	800μl
25mM MgCl ₂	500μl	500μl
6×loading dye	300 μ l	500 μ l
ddH ₂ O	1250μl	1250μl×2

Description

BioReady Taq is a highly pure, thermostable recombinant DNA polymerase encoded by a modified gene from *Thermus aquaticus* species and expressed in E.coli. Its recombinant nature ensures utmost purity, reproducibility and processivity. BioReady Taq processes an initial 5'→3' exonuclease activity and lacks 3'→5' exonuclease activity. The enzyme leaves a single 3'-nucleotide overhang that make the products suitable for cloning by TA vector system. BioReady Taq provides a thermostability that meets the requirements of specialized PCR applications.

Physical purity

BioReady Taq is determined to be >90% pure as judged by SDS-PAGE gels with Coomassie brilliant blue staining.

Endonuclease Assay

One microgram of super coiled plasmid DNA (pUC18) is incubated with 5 units of BioReady Taq for 8 hours at 45°C and 8 hours at 70°C in 1× reaction buffer. Following incubation, the DNA is visualized as a single band an ethidium bromide-stained agarose gel.

Exonuclease Assay

One microgram of lambda DNA and 1ug lambda/HindIII DNA are incubated with 5 units of BioReady Taq for 8 hours at 45°C and 8 hours at 70°C in 1× reaction buffer. Following incubation, the DNA is visualized as a single band an ethidium bromide-stained agarose gel.

Functional Assay

BioReady Taq is tested for performance in the polymerase chain reaction using 2.5 units of enzyme to amplify a lambda DNA (0.5kb, 1.0kb and 3kb) and beta-globin partial fragment (408bp) of human. The resulting PCR product is visualized as a single band on ethidium bromide –stained agarose gel.

Unit Definition

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTPs into acid-insoluble form within 30 min at 72°C TAPS assay condition.

Storage buffer

20 mM Tris-HCl (pH8.0 at 25°C), 100 mM KCl, 0.1mM DDT, 50% Glycerol, 0.5% Tween 20 and 0.5%NP-40.

Reaction buffer: 10X Reaction buffer with magnesium [proprietary formulation]

Applications

Routine PCR, TA cloning, RT-PCR, Multiplex PCR, PCR-SSCP, PCR-RFLP and PCR-RAPD, etc.

Storage temperature

Store at -20°C. Avoid exposure to frequent temperature changes.

General reaction mix by combining the indicated components:

10 x Buffer	10 µl
MgCl ₂ (25mM)	X µl
dNTP Mixture	8 µl
Sense primer (10 µM)	1-5 µl
Antisense primer (10 µM)	1-5 µl
BioReady Taq	0.5 µl (2.5 units)
Template DNA	10-500 ng
ddH ₂ O	up to 100 µl

General reaction condition:

94°C	3-5min	} 30-40cycles
94°C	30s	
50-65°C	30s	
72°C	1 kb / min	
72°C	5 min	

Company Information

Manufacturer: Hangzhou Bioer Technology Co.,Ltd

Address: No.1192 Bin'An Rd, Binjiang District, Hangzhou, Zhejiang Province,China

Tel: 0571-87774567 Fax: 0571-87774553

Web: www.bioer.com.cn

Zip Code: 310053

Aftersales Service Provider: Hangzhou Bioer Technology Co.,Ltd